


**PLASTIC SURGERY
INSTRUMENTS**



*SURGICAL INSTRUMENTS
MANUFACTURING COMPANY*

DELIVERING PRODUCTS TO THE HEALING HANDS


RETRACTORS	-----	3-16
FIBER OPTIC RETRACTORS	-----	3-10
BREAST DISSECTORS	-----	17-18
MARKERS & CALIPERS	-----	19-20
SPECULAS & SEPTUM FORCEPS	-----	21
ENDO DISSECTORS	-----	22-29
FIBER OPTIC LIGHT GUIDE CABLE	-----	30
PRODUCT INDEX	-----	31



TEBBETS Style Retractors

With Serrated Tips

With Light Guide	With Light Guide & Suction	Blades
SR70-250L	SR70-254SL	16 mm x 8 cm
SR70-251L	SR70-255SL	24 mm x 9 cm
SR70-252L	SR70-256SL	24 mm x 10 cm
SR70-253L	SR70-257SL	25 mm x 10 cm
SR70-258L	SR70-262SL	25 mm x 12 cm
SR70-259L	SR70-263SL	23 mm x 15 cm
SR70-260L	SR70-264SL	30 mm x 9 cm
SR70-261L	SR70-265SL	30 mm x 12 cm
SR70-266L	SR70-269SL	30 mm x 15 cm
SR70-267L	SR70-270SL	36 mm x 16 cm
SR70-268L	SR70-271SL	36 mm x 19 cm



FERRIERA Style Retractors

With Smooth Tips

With Light Guide	With Light Guide & Suction	Blades
SR70-272L	SR70-276SL	14 mm x 7 cm
SR70-273L	SR70-277SL	14 mm x 9 cm
SR70-274L	SR70-278SL	16 mm x 8 cm
SR70-275L	SR70-279SL	16 mm x 15 cm
SR70-280L	SR70-285SL	23 mm x 10 cm
SR70-281L	SR70-286SL	23 mm x 15 cm
SR70-282L	SR70-287SL	23 mm x 18 cm
SR70-283L	SR70-288SL	25 mm x 10 cm
SR70-284L	SR70-289SL	25 mm x 12 cm
SR70-290L	SR70-294SL	25 mm x 15 cm
SR70-291L	SR70-295SL	25 mm x 19 cm
SR70-292L	SR70-296SL	30 mm x 12 cm
SR70-293L	SR70-297SL	30 mm x 15 cm




BIGGS Mammoplasty Retractors

Retractor Only	With Light Guide	With Light Guide & Suction	Blades
SR70-298	SR70-300L	SR70-302SL	30 mm Narrow
SR70-299	SR70-301L	SR70-303SL	50 mm Wide



CIRCULAR BREAST RETRACTOR

SR70-309
WITHOUT FIBER OPTIC
OR SUCTION



CIRCULAR BREAST RETRACTOR

SR70-310
WITH SUCTION
2" arm with
recessed suction tube



CIRCULAR BREAST RETRACTOR

SR70-311
WITH FIBER OPTIC
AND SUCTION



**With Endoscopic Channel & Stopcock,
to hold 10 mm Telescopes.**

SR70-312 25 mm x 18 cm



Mutli Approach Retractor
With Smooth Round Tip

Blades

SR70-313 23 mm x 18 cm



Face-Lift Retractor

Retractor Only	With Light Guide	Blades
SR70-314	SR70-307L	22 mm x 11 cm

WIEDER Nasal Retractor

With Light Guide	With Light Guide & Suction	Blades
SR70-328L SR70-329L	SR70-330SL SR70-331SL	38 mm Small 45 mm Large



WIEDER



MALINIAC

MALINIAC Nasal Retractor

Length 4" (10 cm)

With Light Guide	Blades
SR70-315L	9 mm x 45 mm



CONVERSE

CONVERSE Nasal Retractor

Length 3½" (9 cm)

Retractor Only	With Light Guide	Blades
SR70-316	SR70-320L	10 x 20 mm
SR70-317	SR70-321L	12 x 30 mm
SR70-318	SR70-322L	12 x 40 mm
SR70-319	SR70-323L	12 x 50 mm



AUFRICHT


AUFRICHT Nasal Retractor

Retractor Only	With Light Guide	Blades
SR70-324 SR70-325	SR70-326L SR70-327L	6¼" (16 cm) x 7 mm 7" (18 cm) x 10 mm

HARRINGTON Retractors

12½" (32 cm)

Retractor Only	With Light Guide	With Light Guide & Suction	Blades
SR70-355	SR70-357L	SR70-359SL	40 x 127 mm
SR70-356	SR70-358L	SR70-360SL	60 x 127 mm



HARRINGTON

DEAVER



FREEMAN

FREEMAN Flap Retractors

Retractor Only	With Light Guide	Blades	Length
SR70-332	SR70-335L	13 mm	7" 18 cm
SR70-333	SR70-336L	12 mm	8" 20 cm
SR70-334	SR70-337L	25 mm	8" 20 cm

DEAVER Retractors

Retractor Only	With Light Guide	Blades	Length
SR70-338	SR70-340L	3/4" 19 mm	7" 18 cm
	SR70-341L	1/2" 13 mm	8" 20 cm
SR70-339	SR70-342L	7/8" 22 mm	8" 20 cm
	SR70-343L	1" 25 mm	8" 20 cm
SR70-344	SR70-348L	1" 25 mm	9" 23 cm
SR70-345	SR70-349L	1" 25 mm	12" 30 cm
SR70-346		1" 25 mm	13" 33 cm
SR70-347		1" 25 mm	14" 35 cm
SR70-350	SR70-353L	1½" 38 mm	12" 30 cm
SR70-351		2" 50 mm	12" 30 cm
SR70-352	SR70-354L	3" 75 mm	12" 30 cm



SIGMOID NOTCH



SIGMOID NOTCH
Extra Oral Retractors
Length 8¾" (22 cm)

Retractor Only	With Light Guide	Working Tip
SR70-361	SR70-362L	9.5 / 16 mm



PETRI PTERYGOID




PETRI PTERYGOID
Intra Oral Retractors
Length 9" (23 cm)

Retractor Only	With Light Guide	Working Tip
SR70-363	SR70-364L	25.5 mm

FOR RIGHT CHEEK

FOR LEFT CHEEK



BAUER
Intra Oral Retractors
Length 8¼" (21 cm)

Retractor Only	Working Tip
SR70-365	26.5 cm / 16 mm
SR70-366	26.5 cm / 16 mm

FOR RIGHT CHEEK

FOR LEFT CHEEK



BAUER
Intra Oral Retractors
Length 8¼" (21 cm)

With Light Guide	Working Tip
SR70-367L	26.5 cm / 16 mm
SR70-368L	26.5 cm / 16 mm



LEVASSEUR MERRIL



SHEA

LEVASSEUR MERRIL

Retractors For Vertical Osteotomies
Length 9¾" (25 cm)

Retractor Only
SR70-369

With Light Guide
SR70-370L

Working Tip
12.5 / 27 mm

SHEA

Multipurpose Retractors
Length 9¾" (25 cm)

With Light Guide
SR70-371L

Working Tip
14.5 mm



BREAST SPATULA
With Handle
Length 12½" (31.5 cm)

Working Tip

SR70-372	25 mm x 15.5 cm
SR70-373	38 mm x 15.5 cm



BREAST RETRACTOR
Double-Ended
Length 12½" (31.5 cm)

Working Tip

SR70-374	22 mm x 55 mm
	38 mm x 70 mm

FREEMAN Free-Lift Retractor

With Bent Neck
Length 7" (18 cm)

Blade
SR70-375 38 mm



FREEMAN
Free-Lift Retractor



FREEMAN
Rake Retractor



22 mm



22 mm
2 Recessed Prongs



38 mm



38 mm
2 Recessed Prongs

FREEMAN Rake Retractor

Length 7" (18 cm)

	Blades	Prongs
SR70-376	22 mm	4
SR70-377	22 mm	4 (2 Recessed)
SR70-378	38 mm	4
SR70-379	38 mm	4 (2 Recessed)



FREEMAN
Rake Retractor

FREEMAN Rake Retractor

Length 4" (10 cm)

Blade
SR70-380 22 mm
SR70-381 38 mm



GORNEY Rake
Retractor

GORNEY Retractor

With Spring Neck
Length 3" (7.5 cm)

Blade
SR70-382 38 mm




MAXWELL Flap
Retractor

MAXWELL Flap Retractor

Length 4" (10 cm)

Blade
SR70-383 22 mm
SR70-384 38 mm



CORONAL Lift Retractor

Length 7" (18 cm), 7 Prongs

	Blade	Grip Hook
SR70-385	70 mm	Front
SR70-386	70 mm	Reverse

Reverse Grip

Front Grip

Coronal Brow



CORONAL Lift Retractor
Length 8¾" (21 cm), 6 Prongs


SR70-387 **Blade Width**
6" 15 cm



YANCOSKIE Abdominoplasty Retractor
Length 7½" (19 cm), 6 Prongs

SR70-388 **Blade Width**
6" 15 cm

MILLARD



Millard Thimble Hooks

	Blade	Hook
SR70-530	18 mm	Single Hook
SR70-531	18 mm	Double Hook

FREEMAN

Double Skin Hook, Sharp

Length

SR70-536 4" 10cm




FREEMAN



ANDERSON

ANDERSON Bear Claw Retractor

SR70-532	SR70-534	2¾" 7cm	3 Sharp Prongs
SR70-533	SR70-535	4½" 11cm	5 Sharp Prongs



MARTINEZ KILNER

MARTINEZ KILNER Retractors

2 Sharp Hooks, With Thimble

	Length	Tip
SR70-538	¾" 8 cm	10 mm (Swivel)
SR70-539	½" 11.5 cm	7.5 mm

CONVERSE Alar Retractor

Double Ended
Length 4¾" (12cm)

Blade

SR70-537 10 mm / 14 mm Blunt



CONVERSE



CLODIUS Orbital Retractor
Double Ended

Length
SR70-399 13 mm / 20 mm



WISE Orbital Retractor
Double Ended

Length
SR70-400 6½" 16.5 cm



SHEEN

SHEEN Nasal Retractor

Length
SR70-401 6¼" 16 cm



SHEEN

SHEEN Retractor

Length
SR70-402 6" 15 cm



BERNSTEIN

BERNSTEIN Retractor
Double Ended

Blade
SR70-403 8 x 46 mm



CAROLIN
Finger Retractor
Double Ended

Length
SR70-404 4¼" 11 cm



CRONIN
Lip & Nostrill Retractor
Double Ended

Length
SR70-405 6¼" 15.5 cm



CRONIN
Cheek Retractor
Double Ended

Length
SR70-406 6½" 16.5 cm



FISCHL
shoe Horn Retractor
Double Ended

Length
SR70-407 6¼" 16 cm

Spatulated Dissector

Length
SR70-408 13" 33 cm



RIGHT


LEFT

AGRIS - DINGMAN Breast Dissector

Graduated, Length 14¼" 36 cm

SR70-409
SR70-410
SR70-411

PAIR
LEFT
RIGHT



REYNOLDS Transauxillary Breast Dissector

Fenestrated

SR70-412
SR70-413

Length
13" 33 cm
16" 40 cm



McCOLLUM-DINGMAN Submammary Dissector

Length
SR70-414 12" 30 cm




McCOLLUM-DINGMAN Breast Dissector

Length
SR70-415 16¼" 41 cm



Angulated Blades



Breast Dissector
Angulated Blades

	Length	
SR70-416	13" 33 cm	Rigid
SR70-417	16½" 42 cm	Rigid
SR70-418	13" 33 cm	Malleable
SR70-419	16½" 42 cm	Malleable

Spatulated Blades



Breast Dissector
Spatulated Blades

	Length	
SR70-420	13" 33 cm	Rigid
SR70-421	16½" 42 cm	Rigid
SR70-422	13" 33 cm	Malleable
SR70-423	16½" 42 cm	Malleable

Blunt Blades



Breast Dissector
Curved, Blunt Blade

	Length
SR70-428	12" 30 cm

Breast Dissector
Oval Spatulated Blades

	Length	
SR70-424	13" 33 cm	Rigid
SR70-425	16½" 42 cm	Rigid
SR70-426	13" 33 cm	Malleable
SR70-427	16½" 42 cm	Malleable



PITANGUY
Flap Grasping Forceps

Length
SR70-430 5½" 14 cm



PITANGUY
Facial Demarcator

Length
SR70-431 6¼" 16 cm



D'ASSUMPCAO
Face-Lift Marker
Length 6¼" (16 cm)

Ratchet
SR70-432 Without Ratchet
SR70-433 With Ratchet




McKISSOCK Areola Markers

Size
SR70-438 38 mm
SR70-439 42 mm
SR70-440 45 mm
SR70-441 50 mm

ACKERMANN
Areola Markers

Size
SR70-436 5½" 14 cm
SR70-437 6¼" 16 cm




MARTEN
Flap Marker
Length 5½" (14 cm)

Ratchet
SR70-434 Without Ratchet
SR70-435 With Ratchet



McKISSOCK (COLLYER) Caliper

Size
SR70-442 12" 30 mm




HALPREN Breast Caliper

Size
SR70-443 34 mm
SR70-444 38 mm
SR70-445 42 mm



CASTROVIEJO Caliper
Length 3½" (9 cm)

Range
SR70-446 0-15 mm
SR70-447 0-20 mm



COTTLE Cartilage Crusher

With Clip

SR70-448



KILLIAN Nasal Speculus

Length 5" 13 cm

	Light Guide	Length
SR70-449	SR70-453L	35 mm
SR70-450	SR70-454L	50 mm
SR70-451	SR70-455L	75 mm
SR70-452	SR70-456L	90 mm



ASCH Septum Straightening Forceps

	Size
SR70-457	9" (23 cm)



WALSHAM Septum Straightening Forceps

Length 9" (23 cm)

SR70-458	Straight
SR70-459	Left
SR70-460	Right



SCALP Elevators
Length 9¾" (24 cm)

SCALP Elevators
Length 9¾" (24 cm)

		Blade	Shape
A	SR70-461	12.5 mm Wide	Slightly Curved
B	SR70-462	12.5 mm Wide	"S" Shaped

		Blade	Shape
C	SR70-463	15 mm Wide	Slightly Curved
D	SR70-464	12.5 mm Wide	Slightly Curved



E



F



G



H

SCALP Elevators
Length 9¾" (24 cm)

		Blade	Shape
E	SR70-465	7 mm Wide	Straight
F	SR70-466	7 mm Wide	Slightly Curved
G	SR70-467	7 mm Wide	Strongly Curved

SCALP Elevators with Suction
Length 9¾" (24 cm)

		Blade	Shape
H	SR70-468	7 mm Wide	Slightly Curved
Length 10¼" (26 cm)			
		Blade	Shape
H	SR70-469	8.5 mm Wide	Slightly Curved



FRONTOTEMPORAL Dissectors

Length 9" (22.5 cm)

		Blade	Shape
A	SR70-470	5 mm	Slightly Curved

Length 9" (23 cm)

		Blade	Shape
B	SR70-471	7 mm	Slightly Curved

Length 9" (23 cm)

		Blade	Shape
C	SR70-472	10 mm	Slightly Curved

FRONTOGLABELLAR Dissectors

Length 9" (22.5 cm)

		Blade	Shape
B	SR70-473	7 mm Tip	"S" Shaped



A



B



C



D

Flap Dissectors

Length 9½" (23.5 cm)

		Size	Shape
A	SR70-474	5 mm Tip	Straight

Orbital Rim Dissectors

Length 9" (22.5 cm)

		Size	Shape
C	SR70-476		Slightly Curved

Zygomatic Arch Elevator

Length 8¾" (20.5 cm)

		Size	Shape
B	SR70-475	5 mm Tip	Curved Tip

Temporal T-Dissectors

Length 9½" (23.5 cm)

		Size	Shape
D	SR70-477	10 mm	Slightly Curved



A



B



C

Periosteal Elevators

	Size	Blade	Shape
A	SR70-478 10" 25 cm	12 mm	"S" Shape

Periosteal Spreaders

	Size	Blade	Shape
B	SR70-479 9¾" 24 cm	8 mm Wide	"S" Shape
C	SR70-480 9¾" 24 cm	10 mm Wide	"S" Shape



Arcus Marginalis Dissectors
Length 9¾" (24 cm)

		Blade	Shape
A	SR70-481	10 mm Tip	"S" Shape

Transoral Dissectors
Length 9¾" (23.5 cm)

C	SR70-483
---	----------

Midface Facial Dissectors
Length 9¾" (24 cm)

		Blade	Shape
B	SR70-482	12 mm Tip	Multiple Curves



Corrugated & Procerus Muscle Dissector

Length 7¾" (19.5 cm)

		Blade
A	SR70-484	4.7 mm

Nerve Protector

Length 6¾" (17 cm)

		Blade	Shape
B	SR70-485	4.5 mm	Curved Left
C	SR70-486	4.5 mm	Curved Right

Nerve Cutting Knife

Length 6¾" (17cm)

		Blade	Shape
D	SR70-487	3 mm	Curved

Knife Handle for Plastic Surgery

Length 8¾" (20.5 cm)

E	SR70-488
---	----------

Nerve Hook

Length 9¾" (24 cm)

		Shaft
F	SR70-489	Rigid
	SR70-490	Malleable



A

Nerve Dissector
Length 9¾" (24 cm)

		Blade
A	SR70-491	4.7 mm

Left



B

Nerve Hook
Length 9¾" (24 cm)

		Hook
B	SR70-492	Left
C	SR70-493	Right

Right




C



Light Guide Cables (Without Adaptors)

Black	Gray			
SR70-494	SR70-498	6	feet 183 cm	ø 3.5 mm
SR70-495	SR70-499	7.5	feet 229 cm	ø 3.5 mm
SR70-496	SR70-500	6	feet 183 cm	ø 5.0 mm
SR70-497	SR70-501	7.5	feet 229 cm	ø 5.0 mm



Light source end adaptors for cables

Adaptor	Types
SR70-502	Olympus
SR70-503	Storz
SR70-504	ACMI British Standard
SR70-505	Wolf

Light source end adaptors for cables

Adaptor	Types
SR70-506	Olympus
SR70-507	Storz
SR70-508	ACMI British Standard
SR70-509	Wolf

A	
Adaptors for Cables	30
Instruments End Adaptor	30
Light source end adaptor	30
Arcus Marginalis Dissectors	30
B	
BAUER Retractor	09
BERNSTEIN Retractor	15
BIGGS Mammoplasty Retractors	03
Breast spatula	11
BREAST Retractor (Double-Ended)	11
BREAST Dissectors	17-18
C	
Calipers	20
CAROLIN Finger retractor	15
CASTROVIEJO Caliper	20
CLODIUS Orbital Retractor	15
CONVERSE Nasal retractor	06
CORONAL Brow Lift retractor	12
CORONAL Lift Retractor	13
CONVERSE Alar Retractor	14
COTTLE Cartilage Crusher	21
Corrugated & Procerus muscle dissector	28
CRONIN Lip & Nostrill retractor	16
CRONIN Cheek retractor	16
D	
D'ASSUMPCAO Face-lift marker	19
DEAVER retractor	07
E	
ENDO DISSECTORS	22-29
Arcus Marginalis Dissectors	27
Corrugated & Procerus muscle dissector	28
FRONTOTEMPORAL dissector	24
FRONTOGLABELLAR Dissectors	24
Flap Dissectors	25
Knife handle for Plastic surgery	28
Midface Facial Dissectors	25
Nerve Protector	28
Nerve cutting knife	28
Nerve dissector	29
Nerve Hook	28-29
Orbital Rim Dissectors	25
Periosteal Elevators	26
Periosteal Spreaders	26
SCALP Elevators	22-23
Temporal T-dissectors	26
Transoral dissectors	27
Zygomatic arch elevator	25
F	
Face-lift Retractor	22-29
FERRIERA Style Retractor	03
Fiber optic Light Guide	30
Light Guide cables	30
Light Source End Adaptors	30
Instruments End adaptor	30
FISCHShoe Horn Retractor Double Ended	22-29
Flap dissector	25
FREEMAN Flap Retractor	07
FREEMAN Face-Lift Retractor	12
FREEMAN Rake Retractor	12
FREEMAN Hook retractor	14
FREEMAN Areola markers	19
Frontotemporal Dissectors	24
Frontoglabellar	24
G	
GORNEY Rake Retractor	12

H	
HARRINGTON Retractors	07
HALPREN Breast Caliper	20
K	
KILLIAN Nasal Speculas	21
KNIFE Handle for Plastic Surgery	28
L	
LEVASSEUR Merrill retractor	10
M	
MALINIAC nasal Retractor	06
Markers	19
MARTINEZ Kilner Retractor	14
MARTEN Flap Marker	19
MAXWELL Flap Retractor	12
MCCOLLUM-Dingman Submammary Dissector	17
MCCOLLUM-Dingman Breast Dissector	17
MCKISSOCK areola Markers	19
MCKISSOCK (Collyer) Caliper	21
Midface Facial Dissectors	27
Millard Thimble hooks	14
Multi approach Retractor	05
N	
Nasal speculas	21
Nerve Protector	28
Nerve cutting Knife	28
Nerve hook	28-29
Nerve Dissector	29
O	
ORBITAL Rim Dissectors	25
P	
PETRI PTERYGOID Retractor	08
Priosteal Elevators	26
Priosteal Spreaders	26
PITANGUY Flap Grasping Forceps	19
PITANGUY Facial demarcator	19
R	
Retractor	4-16
Fiber Optic Retractors	4-10
Intra Oral Retractors	8-9
Extra Oral Retractors	08
REYNOLDS Transauxillary Breast dissector	17
S	
Scalp Elevators	22-23
SIGMOID NOTCH retractors	08
SHEA retractors	10
SHEEN Nasal retractor	15
SHEEN Retractors	15
T	
TEBBETS Style Retractors	04
Temporal T-Dissector	25
Transoral Dissectors	27
W	
WALSHAM Septum Straightening Forceps	21
wieder tongue depressors	06
WISE Orbital Retractor	15
Y	
YANCOSKIE Abdominoplasty Retractor	13
Z	
ZYGOMATIC Arch Elevator	25



DELIVERING PRODUCTS TO THE HEALING HANDS

Tax ID: 0105561085653



info@surgyx.com | www.surgyx.com
SurGyx Surgical Company Limited
233/90 Sanpawut Town Sanpawut Road,
Bangna Bangkok, Thailand 10260